


INSOLE-TECHNOLOGY FOR CRAFTSMANSHIP



OUR INSOLE NOVELTIES 2017

INSOLE TECHNOLOGY FOR SPECIAL AREAS OF APPLICATION

SPANNRIT has been standing for innovative solutions in the field of insole supplies for many years now. This year we once again present you the newest insole models in orthopaedic shoe technology, as we have made it our goal to manufacture the right orthopaedic insole for every shoe and every foot ailment and to offer you the respective blanks.

Thanks to the state of the art processing and the selected materials our insoles are also suitable for special areas of application. The modern production technology allows swift order processing, a flexible subsequent blank treatment at your facility and warrants continuous high quality.

MANY BENEFITS

- Good abrasability and bonding properties
- 96% resetting force
- Pre-shaped – from top and bottom
- High moisture absorption
- Optimum pressure distribution
- Machine washable

BALLERINA INSOLE

Size 35 – 44 • available ex stock

Art.-No. M3038B29E3098



- Material:** Open-celled high-grade black MEMOPUR® foam with approx. 18 Shore, padding embedded in foam, upper layer Sahara Skin black, unperforated.
- Composition:** Stability core fully embedded in foam with flexible heel area. The slim build, with heel and fore foot phasing out toward 0, is perfect for use in shoes with very limited space e. g. Ballerinas.
- Properties:** Good edge stability, pre-shaped from top and bottom, high moisture absorption, optimum pressure distribution, machine washable up to 40° C.



PUMPS INSOLE

Size 35 – 42 • available ex stock

Art.-No. M3138B29E3198



- Material:** Open-celled high-grade black MEMOPUR® foam with approx. 18 Shore, extra fore foot padding embedded in foam, upper layer Sahara Skin black, unperforated.
- Composition:** Stability core fully embedded in foam with flexible heel area. Additional padding for the metatarsal head and a high Longitudinal arch.
- Properties:** Good edge stability, pre-shaped from top and bottom, high moisture absorption, optimum pressure distribution, machine washable up to 40° C.



CALCANEAL SPUR INSOLE

Size 35 – 48 • available ex stock

Art.-No. M1678B29E1998



Material: Open-celled high-grade red MEMOPUR® foam with approx. 16 Shore, upper layer Sahara Skin black unperforated. The calcaneal spur padding is protected using a special film that transmits the padding effect in an optimum manner.

Composition: Stability core fully embedded in foam. Cavity in the calcaneal spore area and lowered along the plantar fascia, offering an optimum alleviation in the area of the calcaneal spur.

Properties: Good edge stability, pre-shaped from top and bottom, high moisture absorption, optimum pressure distribution, machine washable up to 40° C.



RIGIDUS INSOLE

Size 35 – 48 • available ex stock

Art.-No. M3238B29E3298



Material: Open-celled high-grade black MEMOPUR® foam with approx. 18 Shore, upper layer Sahara Skin black, unperforated.

Composition: Stability core fully embedded in foam with an extended plate that braces the metatarsophalangeal joint as well as soft padding for the heel.

Properties: Good edge stability, pre-shaped from top and bottom, high moisture absorption, optimum pressure distribution, machine washable up to 40° C.



SUPPORTING INSOLE

Size 35 – 48 • available ex stock

Art.-No. M1638B290098



Material: Open-celled high-grade black MEMOPUR® foam with approx. 18 Shore, upper layer Sahara Skin black, unperforated.

Composition: Stability core fully embedded in foam with integrated soft padding for the heel and a pelotte laid in foam.

Properties: Good edge stability, pre-shaped from top and bottom, high moisture absorption, optimum pressure distribution, machine washable up to 40° C.



NOW ONLINE: THE NEW WEB PORTAL FUSSSCHMERZ-RATGEBER.DE

The new guidebook-portal www.fusschmerz-ratgeber.de describes different ailments of the feet, knees, legs and hip and outlines the benefits of orthopaedic insoles. Users can directly access the link to their local specialised dealer and obtain information. **This service is exclusive to SPANNRIT customers!**

Following contents can be accessed via the portal:

- Address and contact details of listed specialised dealers
- Geographic mapping of listed specialised dealer (Google Maps)
- Foot ailment encyclopedia
- Field reports

MANY BENEFITS

- Easy acquisition of new customers
- More website visitors
- **Listing on the web portal**
- Google friendly via links
- Insole expert image



When visiting www.fusschmerz-ratgeber.de

Potential customers local to your area find your specialised company for orthopaedics.

The map shows all foot experts, that have decided they want to be listed on our web portal. Having your address and a link to your website listed on our web portal is only one click away.

Contact us now and request your application form to have your company listed on our portal: +49 (0) 6027 / 4064-0.

SPANNRIT GmbH
Industriestraße 3
63801 Kleinostheim
Germany
Phone: +49 (0) 6027 / 4064-0
Fax: +49 (0) 6027 / 4064-20
e-mail: info@spannrit.de
www.spannrit.de

